

LightRun Plate

Sample Submission Guide for the 555 Sample Submission Galact LightRun Plate Sequencing Service

OPTIMUM PRIMER CONDITIONS

- The optimum primer length is 16-25 bases
- Primer melting temperature (Tm) should be 50-62° C
- The GC content of the primer should be 35-60 %
- Ideally one G or C should be located at the 3' primer end
- The number of 3' Gs or Cs should not exceed 2 Gs or Cs
- If possible, avoid >3 identical bases in a row in the sequence
- Primers must not contain phosphorylation or fluorescent dyes

There are many DropBoxes installed throughout Europe for free sample shipment.

Eurofins Genomics Sequencing GmbH Gottfried-Hagen-Straße 20 51105 Köln, Germanv



The DNA Universe

Plasmid

Purified

Products

accurate results.

DNA

PCR

DNA CONCENTRATION & VOLUME

150-300 bp

1000-3000 bp

300-1000 bp

80-100 ng/µl

2 ng/µl

12 ng/µl

25 ng/µl

Quantify the DNA concentration via agarose gel or a photometer to ensure

5 µl

5 μl 5 μl

5 ul

Sample Preparation & Submission.





1. PREPARE SAMPLES

Take 5 µl purified template DNA with either of the concentrations given in the adjacent table.

Add 5 μ l primer with a concentration of 5 pmol/ μ l. Ensure that the total volume of your premixed sample is 10 μ l.

2. LABEL PLATES

Use our LightRun 96 Barcodes to label your 96 well PCR plate with V-shaped wells.

Use our LightRun 384 Echo Plate for submitting your samples in 384 well

format.

LightRun 96 Barcodes

LightRun 384 Echo Plate

3. SUBMIT SAMPLES

Seal your 96well plates using 8-cap strips to prevent material loss.

Use our adhesive plate foil to seal the LightRun 384 Echo Plate.

Drop your samples in a Drop-Box to benefit from free sample pick up or ship samples at ambient temperature directly to our sequencing lab.

4. DOWNLOAD RESULTS

You will receive an email once your data is ready to download from your Eurofins Genomics account. Data willbe available the next business day after sample receipt.

Additional Service & Information

If you need the perfect sequencing primer, visit us here! To find your closest DropBox, login to your account and click on "My DropBoxes" Do you need boxes for safe shipment of your sample plates? Visit us here

USE OUR FREE SEQUENCING PRIMER DESIGN TOOL TO DESIGN THE OPTIMUM SEQUENCING PRIMERS

If you search for the perfect sequencing primer, you will love our SeqPrimer and Standard Primer (incl. NightXpress shipping options). Test them with a 25,- Euro voucher (or equivalent national currency). Simply use the promotion code "OLIGOPLANET" during check-out.



