

# Whole Plasmid Sequencing




### **SAMPLE CONCENTRATION & VOLUME**

Plasmid	Regular	2.5 - 25 kb	30 ng/µl	20 μΙ
	Large	25 - 125 kb	50 ng/μl	30 μΙ
	XL	125 - 300 kb	50 ng/μl	50 μΙ

BY DROPBOX There are many DropBoxes installed throughout Europe for free sample shipment.	BY POST: Eurofins Genomics Sequencing Lab Gottfried-Hagen-Straße 20 51105 Köln			

- » Please provide your samples either in nuclease-free water or elution buffer (10 mM Tris, pH 8.5)
- >> Please do not use Nanodrop for your concentration measurement. Spectrophotometric techniques are unreliable for the quantification measurement. Please use a fluormetric method, like Qubit.



# Sample Preparation & Submission.

1. www







### 1. DEFINE REACTIONS

Login to your Eurofins Genomics account. Choose "Whole Plasmid Sequencing" in the order menu.

Define your reactions on the order page and allocate the barcodes to your reactions.

### 2. LABEL SAMPLES

Affix the Whole Plasmid Sequencing barcode on your sample tube / plate. Please don't use parafilm to seal your tubes!

### 3. SUBMIT SAMPLES

Use either one of your nearest DropBox for free sample pick up or a courier of your choice.

Make sure to pack your samples safely.


### 4. DOWNLOAD RESULTS

You will receive an email once your data is ready to download.

Data will be available in your online account.

## **Additional Service & Information**

- To find your closest DropBox, login to your account and click on "my DropBoxes"
- Shipping options for sending samples by mail or courier can be found here
- » Do you need boxes for safe shipment of your samples?
  Visit us here





FOR MORE INFORMATION ON WHOLE PLASMID SEQUENCING APPLICATIONS AND OXFORD NANOPORE TECHNOLOGY, VISIT OUR WEBSITE