

DNA SEQUENCING RESULTS GUIDE

Examining Your Results

This **illustrative guide** briefly explains some of the most frequently encountered sequencing data irregularities in your electropherograms.

A good starting point to correctly evaluate sequencing data is the electropherograms of raw and analysed data files.

Using the **original *.ab1 file** gives you the best chance to precisely compare the raw and analysed data traces.

Several excellent freewares such as Applied Biosystems Sequence Scanner exist to easily examine the electropherograms in your sequence data. For instance, the "**Annotation**" **tab** in the Sequence Scanner software displays signal strength, signal to noise ratio, and even reflects the quality of a sequence in an electropherogram





Successful Sequencing Reactions



Electropherogram of a normal DNA sequencing read

Chromatogram Characteristics

- Well-formed, distinctive single coloured peaks
- Peaks evenly separated
- Absence of background signals

Reasons

- Appropriate template and primer concentrations
- Excellent purity of DNA and template
- Optimum primer design
- Excellent primer quality

Tip!

Use our free sequencing primer design tool and our optimised SeqPrimer to achieve optimum primers.



Failed Sequencing Reactions



Electropherogram of a failed sequencing results

Chromatogram Characteristics

- Absence of clearly defined peaks in raw and analysed data
- Occurrence of excess dye peaks
- Very low signal-to-noise ratios (S/N A/C/G/T: <15)

Reasons

- Insufficient or poor quality template and /or primer
- Primer binding site absent, deleted or mutated

Tip!

Don't despair. Give us a call to get some expert advice



Multiple Sequence Signals



Electropherogram of a multiple sequence signals

Chromatogram Characteristics

- Peaks unevenly spaced
- Overlapping of peaks in raw and analysed data
- Presence of artifacts beneath peaks

Reasons

- Contaminated template or primer
- Poor quality template and /or primer
- Multiple priming

Tip!

Use our optimised and validated SeqPrimer to ensure good primer quality .



Low Signal Strength



Electropherogram with low signal strength and signal die out

Chromatogram Characteristics

- Base calls fade off before the end of the read
- Very low peak heights in raw data trace
- Low signal-to-noise ratios (S/N A/C/G/T: <50)

Reasons

- Suboptimal DNA and primer concentrations
- Impure DNA templates
- Primer mutated
- Poor primer design

Tip!

Use our free sequencing primer design tool for optimal primer design.



Very Strong Signals and Pull-Up Peaks



Electropherogram with very strong signals and Pull-up peaks

Chromatogram Characteristics

- Very high peaks in the raw data trace
- High peaks in the analysed data with pull-up secondary peaks and poor base calls
- Very high signal-to-noise ratios S/N A/C/G/T: > 750

Reason

 Excessive template during cycling sequencing

Tip!

Reduce the amount of template



Mixed Sequenced Reactions

Characteristic of mixed sequence reactions

- More than one sequence in the analysed trace data
- More than one sequence starting after base 30 to 100 (multiple cloning site)
- Lower raw data peaks

Reasons for mixed sequence reactions

- Mixed plasmid preparations
- Multiple PCR products
- Multiple priming sites
- Multiple primers in reaction mix e.g. due to incomplete removal of PCR primers
- Primer-dimer contamination
- Frame shift mutation
- Primer with n-1 contamination or degrading primer
- Slippage due to homopolymer or repeat regions in the template





Mixed Sequences Starting From One Cloning Site



Electropherogram with mixed sequencing starting at a cloning site

Chromatogram Characteristics

• Overlapping peaks

Reason

• Use of more than one colony for plasmid growth and preparation

Tip!

Streak out single colonies at least on a selected medium



Mixed Sequences Starting From The Beginning



Electropherogram with mixed sequencing starting from the beginning

Chromatogram Characteristics

- Overlying sequences throughout the entire read or up a certain region
- Several sequences in the analysed trace data
- More than one sequence starting after base 30 to 100 (multiple cloning site)

Reasons

- Presence of several templates in a sequencing reaction
- Primer dimer contamination

Tip!

Use our free sequencing primer design tool which also checks for primer dimer.



Sequencing Result Comment – Degraded Primer



Electropherogram from degraded primers

Chromatogram Characteristics

- Mixture of different peaks
- Low raw data peaks

Reason

 Unspecific binding due to a mixture of various primer products with differing lengths

- Isolate new DNA from a pure single colony and re-sequence (re-streak the clone if necessary)
- Optimize PCR conditions if multiple bands appear after checking PCR template on an agarose gel
- Gel purify your PCR product
- Ensure primer has one binding site
- Ensure that only one primer is used for sequencing
- Check for remaining primers after PCR clean-up on an agarose gel
- Re-design primer and/or optimise PCR amplification



Long Homopolymer Regions





Electropherogram with homopolymer regions

Chromatogram Characteristics

 Mixed sequences downstream of long homopolymer T or A regions

Reason

 Slippage of the enzyme: incorporation of 20 T bases or 21 or 22 T in a 20 T region. Effect more pronounced in PCR templates than plasmids

- Sequence the complementary strand
- Try using an anchored primer: e.g. polyT sequencing primer with an A, C or G base at the 3' end of the polyT to anchor the primer at the end of the homopolymer region
- Special sequencing conditions (low annealing temperatures) necessary with anchored primers. Use our TubeSeq service and select PowerRead upgrade.



Top Heavy Sequences



Electropherogram with heavy sequences

Chromatogram Characteristics

- Very high peaks at the beginning followed by a sharp decrease in signal intensity
- Short fragments are generated in excess and preferentially injected into the capillary

Reasons

- Template or primers are depleted in the early rounds of cycle sequencing creating an excess of short fragments
- Short fragments are preferentially injected into the capillary leading to blockage which inhibits the rapid flow of reaction products
- Excessive use of template or primer in the sequencing reaction



Abrupt Signal Loss



Electropherogram with abrupt peak loss

Chromatogram Characteristics

- High peaks that stop abruptly at a certain position
- Base calls may stop before the end of a read

Reasons

- Template has a secondary structure
- High GC content

- Sequencing the template from the opposite direction
- Use a primer at a different position
- Apply 7-deaza-dGTP at the PCR amplification stage to relax the hairpin structure (75% of the dGTP in the PCR)
- Use our TubeSeq service and select PowerRead upgrade



Repetitive Regions



Electropherogram with repetitive regions (in this case GT- repeat)

Chromatogram Characteristics

- Gradually decreasing peak height starting from within the repeat
- Base calls are fading off after the repeat region)

Reason

• Repetitive regions may lead to a dropdown in signal intensities

- Sequencing the template from the opposite direction
- Use a primer at a different position