

SAMPLE SUBMISSION GUIDE

Site-directed mutagenesis

Sending samples according to the requirements below helps us to do our job better and provides you with accurate results.

Sample Data and Formats

- Please use our online forms for SDM service quote requests
- We will calculate the price of the project based on your request and send you a quotation. The quotation summarizes your request and specifies the price of the service.

Sample Requirements

- Please send us all samples at a concentration of 100-200ng/μl
- The total sample volume should be minimum 30μl (in total around 3μg)
- The sample should be sent in double distilled water or TE buffer
- Please send your samples in 1.5ml tubes
- Labelling of the tubes:
 - Sample name as indicated in the request form
 - Concentration of the sample
 - Volume of the sample

Shipping

- Please mention the Order ID
- Please use a padded envelope when sending the vector to ensure that it arrives safely
- Ship the samples liquid at ambient temperature
- Shipping address:

**Gene Synthesis Department
(Project Assistant-Lab)
Eurofins Genomics
Anzinger Str. 7a
D-85560 Ebersberg**

Contact

Please do not hesitate to contact us in case of any questions!

Email: genes-eu@genomics.eurofinseu.com

Phone: +49 (8092) 3379800