

96well FrameStar[®] PCR Plate

Thermal Cycler Compatibility Table

Manufacturer	Standard Cyclers (Model type)	qPCR Cyclers (Model type)
Biometra	Uno T1 Thermocycler T Gradient T Robot	
Bio-Rad	PTC100™ (96well block only) PTC200 (DNA Engine™) PRC220/221 (DNA Dyad™) PTC225 (DNA Tetrad™)	CFX96
Bio-Rad-MJ Research		Opticon™ Opticon2™ Chromo4™
Corbett Research (Qiagen)	Palm Cyclers	
Eppendorf	MasterCycler® Gradient MasterCycler® EP Gradient	Mastercycler™ ep realplex
G-Storm	GS1-GS4 (all models) GSX	
MWG	Primus 96 The Q-Lifecycler	
Peqlab	peqStar 96	
SensoQuest	labcycler	
Stratagene	Gradient Cyclers	
Techne	Multiblock System & MBS Omnigene Omni-E PCR Express / Px2 / PxE PCR Sprint Touchdown	